

advantage New Promotion Sept. 1–Dec. 31, 2017



Like It, Order It, Love It!

New special offers for centrifuges, freezers, and more.  
Cool anniversary savings guaranteed!



Centrifuge 5810 R: renowned quality and reliability in a compact footprint! Get it and love it!



Eppendorf ULT freezers: wide choice of options for high performance, sample safety, and convenience

## You'll Simply Love It!

Do you plan to replace an old instrument? Or do you have funds for additional equipment? So you are busy studying product literature and researching the web. You have discussed with sales reps, colleagues and, of course, your purchaser. You may even consider visiting a tradeshow for a live experience.

We can't promise when it will happen, but eventually, you'll fall in love with an Eppendorf instrument! And whatever Eppendorf product you like, there is nothing easier than adding your new favorite to your shopping cart and placing your order in our online shop, or by getting in touch with us personally for competent expert advice.

## 15 Years Eppendorf Nordic

Take advantage now of cool anniversary savings for refrigerated centrifuges, ultra-low temperature freezers, and more!

- > Centrifuge 5810 R Cell Culture/Plate Package with 15% savings
- > Centrifuge 5430 R Tube Volume Variety Package with 15% savings
- > Up to 20% savings on Innova® U725 and U725-G
- > ... and more great anniversary offers for centrifugation products, tubes, and freezer accessories

# Centrifuge 5810 R in a New Package: Save 15 %!

Centrifuge 5810 R with its renowned quality and reliability offers you a very cost-efficient solution for your medium to high-throughput applications. It combines outstanding versatility and capacity for both tubes and plates with an extraordinary compact footprint. Centrifuge 5810 R is the most compact refrigerated 3 liter centrifuge currently available.

The 5810 R Cell Culture/Plate Centrifuge Package offers a capacity of max. 28 x 50 mL/56 x 15 mL and 2 plate buckets for each 1 deepwell plate/4 microplates.

## NEW! Cell Culture/Plate Centrifuge Package

- > Centrifuge 5810 R, refrigerated
- > Rotor S-4-104, incl. round buckets
- > 2 sets of 2 adapters for 15 mL conical tubes (max. 3,197 x g)
- > 2 sets of 2 adapters for 50 mL conical tubes (max. 3,180 x g)
- > 2 plate buckets (max. 2,568 x g)
- > 500 Eppendorf Conical Tubes 15 mL
- > 500 Eppendorf Conical Tubes 50 mL

*More info on Eppendorf Conical Tubes on page 6!*

**From September 1 until December 31, 2017, the new Cell Culture/Plate Centrifuge Package including Centrifuge 5810 R is available with 15% anniversary savings!**

**Tip! Also see our anniversary offers for additional rotors on page 5!**

## NEW! Centrifuge 5810 R Cell Culture/Plate Package



The new Cell Culture/Plate Centrifuge Package allows mixed loading of plates and conical tubes – for increased versatility!

## Centrifuge 5810 R (refrigerated)

- > Centrifuge lid with soft-touch lid closure
- > Low access height of 29 cm for easy loading and unloading of samples
- > Quiet operation to improve your work environment
- > Temperature range: -9 °C to +40 °C
- > FastTemp function for fast pre-cooling
- > ECO shut-off engages after 8 hours of non-use to reduce energy consumption

\*The Centrifuge 5810 R is an *in vitro* diagnostic device according to Directive 98/79/EC of the European Parliament and the Council dated October 27, 1998.

### Order information

#### Description

Cell Culture/Plate Centrifuge Package

#### Order no.

5811 000.908

### Advantage prices (excl. VAT)

DKK 

70.800,00

EUR 

9.300,00

NOK 

85.900,00

SEK 

87.800,00

## NEW! Centrifuge 5430 R Tube Volume Variety Package



The refrigerated Centrifuge 5430 R is now available in a new tube volume variety package. Choice of rotary knobs or keypad.

### Centrifuge 5430 R (refrigerated)

- > Menu-driven, multi lingual operation with large backlit display
- > 5 program keys for easy access to routine programs
- > Eppendorf QuickLock® system for quickly opening and closing the rotor lid
- > Temperature range: -11 °C to +40 °C
- > Continuous cooling maintains a constant temperature after the run has been completed – your samples stay cool
- > FastTemp pro® function for programmable pre-cooling according to day and time

# Centrifuge 5430 R in a New Package: Save 15 %!

The unique refrigerated cross-over Centrifuge 5430 R combines the best features of a microcentrifuge (small footprint) and multipurpose centrifuge (versatility) in one instrument. In a compact size, it can accommodate rotors for a variety of tubes as well as for microplates, spin columns, and 15 mL and 50 mL conical tubes.

The new Centrifuge 5430 R Tube Volume Variety Package includes:

- > Centrifuge 5430 R, refrigerated
  - > Rotor FA-45-30-11, aerosol-tight, with Eppendorf QuickLock rotor lid; for 30 x 1.5/2.0 mL tubes (max. 20,817 x g)
  - > Rotor F-35-6-30, with rotor lid, incl. adapters for 6 x 15 mL conical tubes (max. 7,540 x g) and 6 x 50 mL conical tubes (max. 7,581 x g)
  - > 500 Eppendorf Conical Tubes 15 mL
  - > 500 Eppendorf Conical Tubes 50 mL
- More info on Eppendorf Conical Tubes on page 6!*

\*The Centrifuge 5430 R is an *in vitro* diagnostic device according to Directive 98/79/EC of the European Parliament and the Council dated October 27, 1998.

**From September 1 until December 31, 2017, the new Centrifuge 5430 R Tube Volume Variety Package is available with 15% anniversary savings!**

**Tip! Also see our anniversary offers for additional rotors on page 5!**

#### Order information

#### Centrifuge 5430 R Tube Volume Variety Package

Rotary knob version

Keypad version

Order no.

5428 000.309

5428 000.953

#### Advantage prices (excl. VAT)

DKK 

49.500,00

EUR 

6.500,00

NOK 

60.000,00

SEK 

61.400,00

Anniversary Rotor Offer: Save 15 %!



Handling hazardous samples requires the use of special rotors, e.g.

① A-4-62, ② FA-45-6-30, and ③ FA-45-24-11-Kit

#### Rotor A-4-62 for Centrifuges 5810/R and 5804/R

- > Swing-bucket rotor; max. speed  $3,220 \times g$  (4,000 rpm)
- > Holds tubes and bottles from 0.2–250 mL as well as plates
- > Aerosol-tight caps are optionally available<sup>1</sup>

#### Rotor FA-45-6-30 for Centrifuges 5810/R and 5804/R

- > 6 × 50 mL conical tube rotor with aerosol-tight Eppendorf QuickLock lid; for tubes from 5–50 mL
- > Max. speed in Centrifuge 5810/R:  $20,130 \times g$  (12,000 rpm)

<sup>1</sup>Aerosol-tightness tested by Public Health England, Porton Down, UK.

# Anniversary Offer: Save 15 % on Special Rotors

We introduced our very first centrifuge in 1964 as an integral part of the Eppendorf Microliter System, and revolutionized life science research practices worldwide. To this day, the Eppendorf name remains synonymous with creative design, innovative technology, and reliable performance that stand the test of time.





A broad rotor program complements our broad range of centrifuges to enable outstanding versatility for all of your applications. Here are a just few examples of special rotors that can facilitate the safe and ergonomic handling of difficult or hazardous samples.

**During the promotional period you can purchase rotors A-4-62, A-4-62 MTP (not shown), FA-45-6-30, and FA-45-24-11-Kit with 15 % anniversary discount.**

#### Rotor FA-45-24-11-Kit for Centrifuges 5430 and 5430 R

- > 24-place rotor with higher edge for safe centrifugation of up to 24 spin columns
- > Aerosol-tight Eppendorf QuickLock lid for safe centrifugation of hazardous samples
- > Max. speed:  $19,090 \times g$  (13,200 rpm)

#### Order information

Description	Order no.	Advantage prices (excl. VAT)			
		DKK 	EUR 	NOK 	SEK 
Rotor A-4-62, with 4 rectangular buckets	5810 709.008	10.455,00	1.377,00	12.750,00	13.005,00
Rotor A-4-62, with 4 MTP buckets	5810 711.002	10.030,00	1.317,50	12.155,00	12.495,00
Rotor FA-45-6-30	5820 715.006	12.920,00	1.700,00	15.640,00	15.980,00
Rotor FA-45-24-11-Kit	5427 752.005	4.649,50	612,00	5.644,00	5.763,00

## Do your conical tubes offer secure protection from sample loss?

Conical tubes 15 mL and 50 mL are used in many laboratory protocols, making them two of the most widely used tube formats in life science laboratories. Given their broad range of applications, the security of the seal is paramount. The specially developed screw cap of the Eppendorf Conical Tubes establishes new standards with respect to handling and sample safety.

### Cap tightness

The typical areas of application for conical tubes include centrifugation, mixing, incubation, transport, and storage of the most diverse samples and solutions. A reliable seal is therefore a basic condition for safe and accurate work.

For these reasons, strong emphasis was placed on this aspect during the development of the screw cap for the Eppendorf Conical Tubes. The exceptional tightness of the seal could be demonstrated through comparative experiments such as storage of ethanol, a volatile liquid (Fig. 1).




Fig. 1: Sample loss after 24 h of horizontal storage in different conical tubes 50 mL at -86 °C. N = 50 Eppendorf Conical Tubes 50 mL and 22 tubes of manufacturers G, Fa, C.

Additional information: [www.eppendorf.com/appnote343](http://www.eppendorf.com/appnote343)



### Tube handling and prevention of contamination

The grooved, multi-surface side contour ensures a non-slip grip while simultaneously allowing secure single-handed opening and closing of the Eppendorf Conical Tubes. Experienced users in particular will handle conical tubes in this fashion for swift and efficient work.

Additional user advantages result from the flat, smooth, and light-colored surface of the screw caps.

Whereas other lids, due to dark coloring schemes and uneven or sometimes concave surfaces, do not offer user-friendly labeling opportunities, the smooth, light gray lids of the Eppendorf Conical Tubes allow legible, high-contrast labeling.

Thanks to their flattened sides, the lids may be placed safely upright on the laboratory bench during sample processing. This minimizes the risk of contamination of the inside surface of the lid as well as the danger of confusing lids during simultaneous handling of multiple samples.



For additional information and the opportunity to order free samples, please visit [www.eppendorf.com/conicals](http://www.eppendorf.com/conicals).

Anniversary MiniSpin® Offer: Save 15 %!

# Anniversary Offer: MiniSpin® with 15 % Off

The Eppendorf approach to product development is about listening to our customers' needs, understanding daily laboratory routines and anticipating future trends. The result is a comprehensive line of microcentrifuges that serve the multiple applications you encounter in your lab.

If you just have to spin a few tubes at a time or if you are looking for a really personal centrifuge that is also fun to work with, take a closer look at the famous MiniSpin from Eppendorf! Thanks to the MiniSpin's extremely compact footprint, you can equip each workstation with a "personal" centrifuge.

\*The MiniSpin is an *in vitro* diagnostic device according to Directive 98/79/EC of the European Parliament and the Council dated October 27, 1998.





MiniSpin®: easy-to-use, powerful entry-level personal centrifuge

## MiniSpin®: your personal microcentrifuge

- > Capacity: 12 × 1.5/2.0 mL tubes
- > Adapters available for tubes from 0.2 mL–0.6 mL
- > Additional rotor option for up to 2 × 8-tube PCR strips
- > Max. speed: 12,100 × *g* (13,400 rpm)
- > User-friendly digital display for time and speed
- > Rotor, inner centrifuge lid, and lid latch are made of metal for maximum operation safety
- > Soft-touch lid closure for ergonomic lid locking
- > Automatic lid opening at the end of the run to prevent sample warming and to allow easy access to samples

**From September 1 until December 31, 2017, the MiniSpin personal microcentrifuge with 12 x 1.5/2.0 mL rotor is available with 15% anniversary discount.**

### Order information

Description	Order no.	Advantage prices (excl. VAT)			
		DKK 	EUR 	NOK 	SEK 
Centrifuge MiniSpin®, non-refrigerated, with Rotor F-45-12-11, 230 V/50–60 Hz	5452 000.018	5.899,00	773,50	7.157,00	7.310,00

Anniversary Offer: Save 15 % on Eppendorf Tubes®!



Trust in the original Eppendorf Safe-Lock Tubes, because your samples deserve nothing but the best.

### Eppendorf Safe-Lock Tubes

- > Hinged Eppendorf Safe-Lock lid prevents unintentional lid opening during incubation and storage
- > Polypropylene provides high resistance to chemicals, mechanical stress, and temperature extremes
- > Precise lid sealing for minimal evaporation rates during long-term storage

# Anniversary Offer: Save 15 % on Eppendorf Tubes®

Eppendorf Safe-Lock Tubes are the result of 50 years of continuous improvement and development. The original Eppendorf Safe-Lock Tubes protect your sample with excellent centrifugation stability, feature perfect sealing properties, and reliably prevent sample evaporation.

They are manufactured without the use of slip agents, plasticizers, and biocides – substances that have been shown to leach from plastic consumables into the sample and negatively affect bioassay results.





**Safer than ever:** Tests at 30,000 × g have shown that vessels can break during centrifugation. Eppendorf Tubes are guaranteed to withstand forces up to 30,000 × g, allowing you to safely centrifuge without sample losses.

**From September 1 until December 31, 2017, selected Eppendorf Safe-Lock Tubes are available with 15% anniversary discount.**

### Eppendorf Biopur® purity grade


All consumables available in this highest possible degree of purity are guaranteed to be sterile and free of pyrogens, RNases, DNases, DNA and ATP as well as PCR inhibitors.

#### Order information

Description	Order no.	Advantage prices (excl. VAT)			
		DKK 	EUR 	NOK 	SEK 
Eppendorf Safe-Lock Tubes, 1.5 mL, Eppendorf Quality™, colorless, 1,000 pcs.	0030 120.086	280,50	36,55	340,00	348,50
Eppendorf Safe-Lock Tubes, 2.0 mL, Eppendorf Quality™, colorless, 1,000 pcs.	0030 120.094	382,50	50,15	459,00	467,50
Eppendorf Safe-Lock Tubes, 1.5 mL, Eppendorf Biopur®, colorless, 100 pcs., individually packed	0030 121.589	187,00	24,65	229,50	238,00
Eppendorf Safe-Lock Tubes, 2.0 mL, Eppendorf Biopur®, colorless, 100 pcs., individually packed	0030 121.597	255,00	33,15	306,00	314,50



Save 15 % on Eppendorf LoBind® Tubes



Eppendorf LoBind Tubes maximize sample recovery by significantly reducing sample binding to the surface.

#### DNA LoBind Tubes are recommended for

- > Preparation or storage of DNA and RNA samples
- > Forensic trace analysis
- > Preparing dilution series in quantitative qPCR
- > Sample preparation in next-generation sequencing
- > Creation of genomic or oligonucleotide libraries

# Anniversary Offer: Save 15 % on LoBind Tubes

#### Eppendorf LoBind® products





- > LoBind material guarantees maximum sample recovery for improved assay results
- > Free of surface coating (e.g. silicone) to minimize the risk of sample interference
- > Lot-certified PCR clean purity grade: free of human DNA, DNase, RNase, and PCR inhibitors
- > Available in tube, microplate, and deepwell plate formats for easy up-scaling
- > Precise lid sealing to minimize evaporation

**During the promotional period selected Eppendorf LoBind Tubes are available with 15 % anniversary discount.**

#### Protein LoBind Tubes are recommended for

- > Preparation or storage of protein, peptide, or antibody samples
- > Enzymatic assays – the hydrophilic surface reduces denaturation when it comes into contact with the tube inner wall
- > Prevention of sample loss during storage of virus stock solutions
- > Storage of cell suspensions

#### Order information

Description	Order no.	Advantage prices (excl. VAT)			
		DKK 	EUR 	NOK 	SEK 
DNA LoBind Tubes, DNA LoBind, 1.5 mL, PCR clean, colorless, 250 tubes (5 bags × 50 tubes)	0030 108.051	85,00	11,05	102,00	102,00
DNA LoBind Tubes, DNA LoBind, 2.0 mL, PCR clean, colorless, 250 tubes (5 bags × 50 tubes)	0030 108.078	102,00	13,60	127,50	127,50
Protein LoBind Tubes, Protein LoBind, 1.5 mL, PCR clean, colorless, 100 tubes (2 bags × 50 tubes)	0030 108.116	59,50	8,50	76,50	76,50
Protein LoBind Tubes, Protein LoBind, 2.0 mL, PCR clean, colorless, 100 tubes (2 bags × 50 tubes)	0030 108.132	102,00	12,75	119,00	127,50

Save 20 % on Eppendorf ULT Freezers



# Eppendorf ULT Freezers: Cool Savings!

Reliability, sustainability, economics, safety – these are the corner stones for selecting the next ULT freezer. And this is exactly what you get from an Eppendorf ULT freezer.



Your new Eppendorf ULT freezer brings high performance, sample safety, and convenience to your lab. You'll simply love it!

**From September 1 until December 31, 2017, Eppendorf ULT freezer models Innova® U725 and Innova U725-G are available with 20 % savings!**





## Innova® U725-G

- > 725 L capacity
- > Holds up to 50,400 samples
- > Floor space needed: 102.5 x 86.7 cm
- > Low energy consumption of 12.2 kWh/day by hydrocarbon-based refrigerants for low running costs
- > Also available as water-cooled version

## Innova® U725

- > 725 L capacity
- > Holds up to 50,400 samples
- > Floor space needed: 102.5 x 86.7 cm

### Order information

Description	Order no.	Advantage prices (excl. VAT)			
		DKK 	EUR 	NOK 	SEK 
Innova® U725	U9440-0001	76.640,00	10.160,00	88.560,00	95.040,00
Innova® U725-G, Water-Cooled	U9440-0004	82.880,00	10.960,00	95.840,00	102.800,00
Innova® U725-G, Air-Cooled	U9440-0005	79.760,00	10.560,00	92.160,00	98.880,00

Take Your Sample Freezing + Storage to the Next Level!



Organized Sample Storage: Our metal freezer racks and polypropylene storage boxes enable optimal use of space.

## Freezers + Accessories: Cool Anniversary Offers

Make your scientific life easier: Store your freezer storage boxes within organized freezer racks. These metal racks optimize your freezer space and make it easy to find the items you desire. Don't waste time sifting through unsorted samples.

Our two cool anniversary offers enable you to take your sample freezing and storage to the next level!





**Offer 1: Tailor-made freezer bundle. Save 15 % on Innova® U725 (see page 10) plus get additional 50 % off a matching Eppendorf Freezer Rack.**

**Offer 2: Get 15 % off selected Eppendorf Storage Boxes.**



Eppendorf Storage Boxes provide safe storage of your sample tubes in a space-saving design.

### Order information

Description	Order no.	Advantage prices (excl. VAT)			
		DKK 	EUR 	NOK 	SEK 
Bundle: Innova® U725 + Freezer Rack for Innova® U725: drawer rack, for 53 mm/2 in boxes, max. 4 x 7 boxes per rack, stainless steel	U9440-0001 + 6001 012.210 <i>(Offer not available in the eShop)</i>	82.085,00	10.880,00	94.850,00	101.790,00
<b>Polypropylene storage boxes, -86 °C, autoclavable, with lid and alphanumeric code; 53 mm/2 in high</b>					
Storage Box 10 × 10, for 100 cryogenic tubes w. int. thread, 3 pcs.	0030 140.508	153,00	20,40	187,00	195,50
Storage Box 9 × 9, for 81 screw cap (cryogenic) tubes 1–2 mL, 3 pcs.	0030 140.516	161,50	21,25	195,50	204,00
Storage Box 8 × 8, for 64 tubes 1–2 mL, 3 pcs.	0030 140.524	161,50	21,25	195,50	204,00

## Expand Your Knowledge

The new online platform Eppendorf Handling Solutions is a comprehensive source of expert knowledge which helps you master the daily challenges of laboratory life.

Discover our Sample Handling solutions on this website and learn more about centrifugation as well as user and sample safety!

[www.eppendorf.com/centrifugationexperts](http://www.eppendorf.com/centrifugationexperts)

**Your local distributor: [www.eppendorf.com/contact](http://www.eppendorf.com/contact)**  
Eppendorf Nordic A/S · 2970 Hørsholm · Denmark  
Tlf. +45 43 24 00 00 · [nordic@eppendorf.dk](mailto:nordic@eppendorf.dk) · [www.eppendorf.dk](http://www.eppendorf.dk)

[www.eppendorf.com/advantage](http://www.eppendorf.com/advantage)

## Takk! Kiitos! Tak! Tack!

Eppendorf Nordic is celebrating its 15<sup>th</sup> anniversary! So, after having your trust over the past 15 years the time has come to thank you.

Celebrate our birthday with us with at least 15% savings on selected Eppendorf products until December 31, 2017.



Eppendorf Norge AS · 0349 Oslo · Norway  
Tlf. +47 22 56 66 32 · [nordic@eppendorf.dk](mailto:nordic@eppendorf.dk) · [www.eppendorf.no](http://www.eppendorf.no)