

Product Description

LoFric Elle, a single use hydrophilic urinary catheter, designed for female users.

Sterilization and shelf life

Sterilized by irradiation. Use before expiry date on the package.

Packaging

30 catheters per customer box. 2 customer boxes per shipping box.

Material composition

Catheter Material:
Polyolefin-based elastomer (POBE).

Surface coating:
Polyvinyl pyrrolidone (PVP).

Wetting fluid:
Water-based wetting fluid with glycerol.

Primary package Material:
Polypropylene.

Product range

| Reference | Type | Length | Size |
|------------|--------------------------------|---------|---------------------------------|
| 418 10 XX* | LoFric® Elle™ Nelaton (female) | 10 cm** | CH 10 |
| 418 12 XX* | LoFric® Elle™ Nelaton (female) | 10 cm** | CH 12 |
| 418 14 XX* | LoFric® Elle™ Nelaton (female) | 10 cm** | CH 14 |
| 418 10 05 | LoFric® Elle™ Nelaton (female) | 10 cm** | CH 10 <small>sample box</small> |
| 418 12 05 | LoFric® Elle™ Nelaton (female) | 10 cm** | CH 12 <small>sample box</small> |
| 418 14 05 | LoFric® Elle™ Nelaton (female) | 10 cm** | CH 14 <small>sample box</small> |


*Generic article number without the two digit suffix specific for a region or country destination when distributing an article.

**Effective length, length of the catheter that can be inserted into the body.

Sample box, 5 pieces.



User-friendly customer box.



- Designed to be ergonomic, user-friendly, with an angulated handle – for good grip and non-touch technique. Connects to a urine collection bag or extension tube.
- Ready to use – Ready when you need it, to enable fast handling and catheterization.
- Slim and Discreet design.
- Lower container is recyclable as plastic (polypropylene).
- Not manufactured with PVC, Phthalates and latex.
- Smooth eyelets.^{1,2}
- LoFric catheters are coated with the unique Urotonic™ Surface Technology. This ensures the catheter remains smooth and safe, both during insertion and withdrawal.^{3,4}

1. Lundgren et al., Spinal Cord, 2000
2. Newman et al., Springer, 2018

3. Waller et al., Spinal Cord, 1997
4. Dunger et al., Acta Endocrinologica, 1988

At Wellspect we develop innovative continence care solutions that change people's lives. We are committed to inspire our users to build self-confidence and independence as well as good health and well-being. We have been leading the industry for over 30 years with our product brands LoFric® and Navina™. We create reliable and user-friendly products for bladder and bowel management with as little environmental impact as possible. We passionately strive to become climate neutral and work closely together with users and healthcare professionals who constantly inspire us to improve our products and services in a sustainable way, now and for the future.

Wellspect. A real difference.

For more information about our products and our initiative Advancing Continence Care Together (ACCT), please visit [Wellspect.com](https://www.wellspect.com).

Join the conversation on Facebook and Instagram.

[wellspect.com](https://www.wellspect.com)

